

# Checklist for Equipment



Ventilator Aspiration with PneuX

## Day 1

### Intubation pack restocking

- 1x size 8 PneuX Tube
- 1x size 9 PneuX Tube
- Stylet
- Inline suction catheter (63cm)
- Instructions for intubation
- Bedside guide

Equipment	
Extension tube	
Cuff pressure monitor	
Inline suction catheter 63cm required?	
Wet circuit (Techs)	
(If delays then 6 hourly saline nebs)	
Bedside Guide	

Check unused PneuX tube size has been returned

## Daily checks

Daily check	
Cuff pressure monitor in use?	
Wet circuit in use?	
Check <u>no</u> HME filter attached if on wet circuit	
Requiring 63cm inline suction?	
Bedside instructions available?	
Check four hourly aspirate on subglottic port	

## Extubation

Cuff pressure monitor returned to store?